

Diagnostic Study of *Rapid Test Chlamydia Ag* in Patients with Genitalia Infection of *Chlamydia trachomatis*

Yuwono

Departement of Microbiology Faculty of Medicine University of
Sriwijaya/Moh.Hoesin General Hospital Palembang South Sumatera Indonesia

Abstract

Genitalia infection of *Chlamydia trachomatis* is a sexual transmitted disease in 10-40% women causing *pelvic inflammatory disease* (PID) with infertility and ectopic pregnancy complications. This infection is very difficult to diagnose due to clinical asymptomatic manifestations. Some of diagnostic methods were developed but still giving unsatisfied results. One of promising method for detecting *C. trachomatis* infection is *rapid test Chlamydia* or *point of care test*. Gold standard of identification *C. trachomatis* was culture or detecting nucleic acid by *polymerase chain reaction* (PCR). This was a diagnostic study conducted in Klinik Graha Sriwijaya Palembang and Departemen Mikrobiologi Klinik RSUP Dr. Moh.Hoesin Palembang from Februari - Mei 2012. Endoservic specimen were collected for Gram staining, *rapid test QuickStripe™ Chlamydia Ag* (Savyon® Diagnostic Israel), DNA extraction and PCR. The 108 subjects with 26.82 years old in average were recruited. Result of QuickStripe™ *Chlamydia Ag* were 23 (21,3%) positive and 85 (78,7%) negative. Result of PCR detection were 21 (19,4%) positive and 87 (80,6%) negative. Value of Kappa test between 2 examiners of *rapid test QuickStripe™ Chlamydia Ag* was 0.893 indicating the examination was reliable. Analysis of *rapid test Chlamydia* result was compared with PCR result showed sensitivity and spesificity 90,4% dan 95,4% respectively (PPV 82,6%; NPV 97,6%; PLR 19,65%; NLR 0,100). The result was indicated that *QuickStripe™ Chlamydia Ag* was good for screening of the disease. Conclusion: *Rapid Test Chlamydia Ag* showed diagnostic significant value and can implementate as a alternative diagnostic method for detecting *C. trachomatis* infection especially in high risk population.

Key Words: *Rapid Test Chlamydia Ag*, PCR, Diagnostic Study

Introduction

The sexual transmitted disease due to *Chlamydia trachomatis* infection is found in 10-40% high risk women.¹ *C. trachomatis* is a intracellular pathogen which destruction of genital mucous layer.² Study by Jazan *et al.* (2003) with subjects women sexual workers in 7 cities in Indonesia showed prevalence 14% - 55% .⁷ Margareth *et al.* (2004) showed that recurrency of the disease up to 40% in 9 months.³ Ascending infection to pelvic area caused *pelvic inflammatory disease* (PID) with infertility and ectopic pregnancy complications.^{2,3,8,9} Some factors contribute to the disease such as age, sexual habit and accuracy of clinical screening. *Centers for Disease Control and Prevention* (CDC) USA recomended screening for women that rare use condom in sexual relation and women with more than 1 partners.⁶

Genital infection of *C. trachomatis* was difficult to diagnose because more than 70% infection was asymptomatic.³ Routin laboratory examination can not different this condition.^{2,12} Some of diagnosis approach have developed but no one satisfaction. Gold standard of diagnosis was culture but it was difficult so identification of nucleic acid by polymerase chain reaction (PCR) can replace the gold standard. PCR is needed a high grade of laboratory. In Indonesia where the laboratory facilities was minial, need another method for diagnosig the infection. A promising method is detection of antibody against lipopolysaccharides of *C. thracomatis* by *rapid test Chlamydia* or *point of care test*. The method was based on chromatography that can give false positive because of cross reaction with other bacteria's lipopolysaccharrides..^{1,15-18} *C. trachomatis* rapid test produced by Savyon[®] Diagnostic (Israel) has 88,5% sensitivity and 96,7% specificity.¹⁸ Tapay *et al.* (2007) got 83,5% and 98,8% in sensitivity and specificity respectively.^{19,1,21} This research is aimed to count diagnostic value of rapid test for detection *C. trachomatis* infection in women sexual worker. It is very important to find promptly diagnosis of asymptomatic manifestation of the disease.

Method

This is an observational study with diagnostic test design. The study conducted in Klinik Graha Sriwijaya Palembang (a clinic in localization of women sexual worker) and Clinical Microbiology Department Moh. Hoesin General Hospital Palembang from Februari until July 2012. The subjects were all women sexual workers in Klinik Graha Sriwijaya Palembang with clinical sign endoservic mucopurulen *discharge*. All the subjects were given informed consent, taking endoservic mucopurulen *discharge* specimen, Gram staining examination, *rapid test Chlamydia* and PCR. Rapid test QuickStripe™ *Chlamydia Ag Savyon® Diagnostic* (Israel) was done according to manufactured guidance. DNA extraction and isolation by *Chelex-100* reaction. PCR using pairs of primer:^{34,53}

Sense: 5'GCCGCTTTGAGTTCTGCTTCC 3' and antisense: 5' GTC GAA AAC AAAGTCACCATA GTA3'. Positive result showed 180 bp in 2% agarose.

Result

The 108 specimen was taken from patients at Klinik Graha Sriwijaya Palembang. Characteristic of subjects were 26.82 years old in average, the youngest 17 years old and the oldest 41 years old. Status of subject's education 60 (55.6 %) elementary school, 32 (29,6%) junior high school and 16 (14.8%) senior high school. Marital status 75 (69,4%) widows, 22 (20,4%) unmarried and 11 (10,2 %) married. Income per month in average Rp.7.330.000 with range Rp.1.000.000,- until Rp. 25.000.000,-. The first sexual intercourse were the youngest 2 (1,9%) in 14 years old and the oldest 2 (1,9%) in 29 years old. Duration of working as women sexual worker were 40 subjects (37,0%) 1 year, 48 subjects (44,4%) 2 years and the rest were 3 – 6 years. Number of customer in last month were 3 persons for 14 subjects (13,0%), 4 persons for 17 subjects (15,7 %), 8 persons for 7 subjects (6,5%), 10 persons for 13 subjects (12,0%) and 15 persons for 12 subjects (11,1%). Compliant using condom during intercourse 83 subjects (76,9%) were rarely, 22 subjects

(20,4.%) routine and 3 subjects (2,8%) never using condom. Clinical complain 62 subjects (57,4 %) symptomatic and 46 subjects (42,5 %) unsymptomatic. Leukorrhoea found in 57 subjects (52,8%) and 5 subjects (4,6 %) disuria.

The result of *rapid test Chlamydia* (QuickStripe™ *Chlamydia* Ag) was shown on table 1.

Table 1. Result of *rapid test Chlamydia*

<i>Rapid test Chlamydia</i>	N	%
Positive	23	21,3
Negative	85	78,7
Total	108	100

Proportion of *C. trachomatis* infection from QuickStripe™ *Chlamydia* Ag was 21,3%.

The result of PCR test was shown in tabel 2

Table 2. Result of PCR test

PCR	N	%
Positive	21	19,4
Negative	87	80,6
Total	108	100

Proportion of *C. trachomatis* infection from PCR test was 19,4%.

Kappa testing from 2 examiners QuickStripe™ *Chlamydia* Ag showed kappa 0.893.

Diagnostic test result were showed in table 3.

Table 3. Diagnostic test results

		PCR		Total
		Positive	Negative	
QuickStripe™	Positive	19 (90,4%)	4 (4,5%)	23(21,2%)
<i>Chlamydia Ag</i>	Negative	2 (9,5%)	83 (95,4%)	85(78,7%)
Total		21(19,4%)	87(80,5%)	108 (100%)

Table 4. Analysis of diagnostic test result

Pemeriksaan	Sn	Sp	PPV	NPV	PLR	NLR
QuickStripe™	90,48	95,40	82,6	97,6	19,68	0,100
<i>Chlamydia Ag</i>						
Sensitivity	= Sn	Spesificity			= Sp	
<i>Positive predictive value</i>	= PPV	<i>Negative predictive value</i>			= NPV	
<i>Positive likelihood ratio</i>	= PLR	<i>Negative likelihood ratio</i>			= NLR	

Discussion

Age average of subject was 26,82 years old with range 17 - 41 years old. Previous research by Krishnaputri in 2011 in the same place age average 25,91 with range 17 - 47 years old.⁵⁴ These data are mostly similar, all subjects in high risk age to get Chlamydia infection. Screening is needed.⁵⁵

Proportion of *C. trachomatis* infection in Palembang by Sedyaningsih et al. (2005) showed 34% higher than this result 19,4%. Although decreasing but still high and need to early diagnostic.⁵⁶ Philips (2006) showed that highly prevalence of *C.*

trachomatis influenced by low grade of socioeconomic and number of costumer as well as in this research.⁵⁵

The first sexuall intercourse lower than 20 years old in this research was a risk factor for *C. trachomatis* infection.⁶ Only 22 subjects (20,4%) routine using condom before intercourse is very low complient as well as study by Sedyaningsih et al. (2005)⁵⁶ and Krishnaputri (2011).⁵⁴ It is also a risk factor for *C. trachomatis* infection.

Most of patients are asymptomatic as well as in more previous research and other research from abroad^{3,54,56}. The data were consistent that Chlamydial infection was asymptomatic and need promptly diagnosis.

Kappa result 0.893 indicate that the *rapid test Chlamydia* is reliable. Analysis from diagnostic test result *rapid test Chlamydia* compare to gold standard PCR showed that number of sensitivity in this result better than result by Saison et al. (2007) in Philippines.²⁰ This high sensitivity suggest that *QuickStripe™ Chlamydia Ag* can be a usefull for screening. Also the specificity result around 96% as well as result by Saison et al. (2007)²⁰. PPV and NPV also support the specificity result. Probability false positive and false negative were also support the specificity as seen in positive and negative *likelihood* ratios. Overall the result from this research showed diagnostic value good compare to gold standard. Thus the rapid test can be used as an alternative method for detecting *C. trachomatis* in high risk population.

DAFTAR PUSTAKA

1. Black CM. Current Methods of laboratory diagnosis of *Chlamydia trachomatis* infections. Clin Microbiol Review 1997; 10: 160-84
2. Stamm WE. *Chlamydia trachomatis* infections of the adult. In: Holmes KK, Mardh PA, Sparling PF, Lemon SM, Stamm WE, eds. Sexually Transmitted Diseases. 3rd New York: McGraw-Hill; 2008: 575-593

3. Margareth R, Hammerschlag. *Chlamydia trachomatis* and *Chlamydia pneumoniae* infection in children and adolescents. *Pediatrics in review* 2004; 25: 43-51
4. Spilopoulou A, Lakiotis V, Vittoraki A, Zavou D, Mauri D. *Chlamydia trachomatis*: time for screening? *Clin Microbiol Infect* 2005; 11: 687-9
5. Honey E, Augood C, Templeton A, Russell I, Pavonen J, Mardh PA, et al. Cost effectiveness of screening for *Chlamydia trachomatis*: A review of published studies, *Sex Transm Infect* 2002; 78: 406-12
6. Peipert JF. Genital chlamydial infections. *N Engl J Med* 2003; 349: 2424-30
7. Jazan S, Tanudya FK, Anartati AS, Gultom M, Purnamawati KA, Sutrisana A, dkk. Laporan Hasil Penelitian Prevalensi Infeksi Saluran Reproduksi pada wanita penaja seks di Bitung Indonesia 2003: 15
8. Bjartling C, Osser S, Jhonson A, Person K. Clinical manifestations and epidemiology of the new genetic variant of *Chlamydia trachomatis*. *Sex Trans Dis* 2009; 36: 529-35
9. Dowe G, Smikle M, King SD, Wynter H, Fredrick J, Hilton-Kong T. High prevalence of genital *Chlamydia trachomatis* infection in women presenting in different clinical settings in Jamaica: implication for control strategies. *Sex Transm Infect* 1999; 75: 412-6.
10. Joesoef MR, Kio D, Linnan M, Kambodji A, Barakbah Y, Idajadi A. Determinant of condom use in female sex worker in Surabaya, Indonesia. *Int J of STD & AIDS* 2000; 11: 262-5
11. World Health Organization. Guideline for management of sexually transmitted infections in female sex workers. Geneva: World Health Organization Regional Office for the Western Pacific 2002: 1-56
12. Crossman, S.H. The challenge of pelvic inflammatory disease. *Am Fam Phsycian* 2006; 73: 859-64
13. Holmes KK, Stamm WE. Lower genital tract infection syndromes in women. In: Holmes KK, Mardh P, Sparing PF, Lemon SM, Stamm WE, Piot P,

Wasserheit JN eds. Sexually Transmitted Diseases 3rd. New York: McGraw-Hill; 2008: 575-93

14. Wahyudi DT, Puspongoro EH, Daili SF. Sensitivitas dan spesifisitas pendekatan sindrom dan jumlah sel polimorfonuklear (PMN) pada infeksi *Chlamydia trachomatis* genital wanita dibandingkan dengan hasil pemeriksaan *Polymerase Chain Reaction* (PCR), 2006. Available at URL: <http://indonesia.digitaljournal.org/index.php/deridn/article/view/26> Diakses tanggal 14 januari 2012
15. Clerc O, Greub B. Routine use of point-of-care test: usefulness and application in clinical microbiology. *Eur J Clin Microbiol Infect Dis* 2010; 16: 1054-61
16. Jensen JS, Bjonelius E, Dohn B, Lidbrink P. Comparison of first void urine and urogenital swab specimens for detection of *Mycoplasma genitalium* and *Chlamydia trachomatis* by polymerase chain reaction in patients attending a sexually transmitted disease clinic. *Sex Trans Dis* 2004; 31: 499-507
17. Rani R, Corbitt G, Killough R, Curless E. Is there any role for rapid test for *Chlamydia trachomatis*? *Int J of STD & AIDS* 2002; 13: 22-4
18. Savyon Diagnostic Ltd. Instruction manual QuickStripe™ *Chlamydia Ag.*
19. Tapay LM, Laitila V, Wawrzyniak, Alexander S, Swain A, Barber P, et al. New point of care *Chlamydia* Rapid test-bridging the gap between diagnostic and treatment: performance evaluation study. *BMJ online first* 2007; 1-8
20. Saison A, Tapay LM, Michel CE, Butress ND, Nadala EC, Magbanua JP et al. Prevalence of *Chlamydia trachomatis* infection among low-and high risk Filipino Women and performance of chlamydia rapid test in Resource-Limited Settings. *J Clin Microbiol* 2007; 45: 4011-17
21. Widjaja S, Cohen S, Brady WE, O'Reilly K, Susanto, Wibowo A, et al. Evaluation of a rapid assay for detection of *Chlamydia trachomatis* infections in out patient clinics in South Kalimantan, Indonesia. *J Clin Microbiol* 1999; 37: 4183-5

22. Van der Pol B, Williams JA, Jones RB. Rapid antigen detection assay for identification of *Chlamydia trachomatis* infection. J Clin Microbiol 1995; 33: 1920-1
23. Martin DH, Nsuami M, Schachter J, Hook III EW, Ferrero D, Quinn TC, et al. Use of multiple nucleic acid amplification test to define the infect patient "gold standard" in clinical trial of new diagnostic test for *Chlamydia trachomatis* infections. J Clin Microbiol 2004; 42: 4749-58
24. Kapala J, Copes D, Sproston A, Patel J, Jang D, Petrich A, et al. Pooling cervical swabs and testing by ligase chain reaction are accurate and cost saving strategies for diagnosis of *Chlamydia trachomatis*. J Clin Microbiol 2000; 38: 2480-3
25. Fadiman KH, Goldman SM. American college of preventive medicine practice policy statement: screening for *Chlamydia trachomatis*. Am J. Prev Med 2003; 24: 287-92
26. Livengood CH, Wrenn JW. Evaluation of Cobas Amplicor (Roche): accuracy in detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* by coamplification of endocervical specimens. J Clin Microbiol 2001; 39: 2928-32
27. Chernesky MA. *Chlamydia trachomatis* diagnostics. Sex Transm Infect 2002; 78: 232-4
28. Vickerman P, Watts C, Allary M, Mabey D, Peeling RW. Sensitivity and *Neisseria gonorrhoeae* in women. Sex Transm Inf 2003; 79: 363-8
29. Schachter J, Stephens RS. Biology of *Chlamydia trachomatis*. In: Holmes KK, Mardh PA, Sparling PF, Lemon SM, Stamm WE, eds. Sexually Transmitted Diseases. 3rd New York: McGraw-Hill; 2008: 555-74

30. Garcia LA, Madkan VK, Tyring SK. Gonorrhea and other venereal diseases. In: Wolff K, Goldsmith LA, Katz SI, eds. Fitzpatrick's Dermatology in General Medicine. 11th. New York; Mc-Graw Hill; 2008: 1996-98
31. Darville T, Hiltke TJ. Pathogenesis of genital tract disease due to chlamydia. J Int Dermatol 2010; 201: 114-25
32. Nyari T, Deak J, Nagy E, Vereb I, Kovacs L, Meszaros G, et al. Epidemiological study of *Chlamydia trachomatis* infection in pregnant women in Hungary. Sex Transm Infect 1998; 74: 213-5
33. Myziuk L, Romanowski B, Brown M, Endocervical Gram stain smears and their usefulness in the diagnosis of *Chlamydia trachomatis*. Sex Transm Infect 2001; 77: 103-6
34. Agha SA, Mashad NE, Rakha SA, Metwally AE. Value of direct antigen detection methods in diagnosing women with suspected genital *Chlamydia trachomatis* infection. African J Microbiol research 2011; 5: 1215-19
35. Sellors JW, Walter SD, Howard M. A new visual indicator of Chlamydial cervicitis? Sex Transm Infect 2000; 76: 46-8
36. Miller KE. Diagnosis and Treatment of *Chlamydia trachomatis* infections. Am Fam Physician 2006; 73: 1411-6
37. Dyck EV, Maheus AZ, Poit P. *Chlamydia trachomatis* infections. In: Laboratory Diagnosis of Sexually Transmitted Diseases. Geneva: WHO Library Cataloguing in Publication Data; 1999: 22-34
38. Peeling RW, Mabey D, Alan Herring and Hook EW. Evaluating diagnostics. Nature Review Microbiology 2006: S7-S19
39. Todd J, Neuwald, Penny L, Urdea M. Biomarker for infection disease diagnostic in the developing world: Diagnosis of *Chlamydia trachomatis* and Neisseria gonorrhoea infection in high risk women. Halteres Associates LLC 2006: 1-28
40. Handojo I Penghantar Imunoasai Dasar. Surabaya: Pusat Penertiban dan Percetakan Unair; 2003: 181-4

41. Schachter J, McCormac WM, Chernesky MA, Martin DH, Van der Pol B, Ricer PA, et al. Vaginal swab are appropriate specimens for diagnosis of genital tract infection with *Chlamydia trachomatis*. J Clin Microbiol 2003; 41: 3784-9
42. Ghanem KG, Jhonson RE, Kauman EH, Marrazo JM, Markowitz LE. Cervical spesimen order and performance measure of *Chlamydia trachomatis* diagnostics testing. J Clin Microbiol 2005; 43: 5295-7
43. Lister NA, Tatrizi SN, Fairley CK, Garland S. Validation of Roche Cobas amplicor Assay for detection of *Chlamydia trachomatis* in rectal and pharyngal spesimen by an amp 1 PCR assay. J Clin Microbiol 2004; 42: 239-41
44. Nester EW, Robert CE, Pearsall NN, Anderson DG. Polymerase Chain Reaction. In Microbiology a Human Perspective, 6th ed. Boston: Mc Graw-Hill 2009: 223-8
45. Jones NL, PCR principles, procedures, and parameters. In: Theopillus BD, Rapley R, eds. PCR Mutation Detection Protocols. Totowa, New Jersey. 2002: 37-46
46. CDC guidelines. Morbidity and Mortality Weekly Report (MMWR). Sex Trans Dis Treatment Guidelines 2010; 55: 38-42
47. Chambers HF. Antimicrobial agent; Protein synthesis inhibitor and miscellaneous antibacterial agents. In: Goodman&Gillman's. The Pharmacological Basis of Therapeutics.10th. New York. McGraw Hill; 2001: 1239-53.
48. Puspinegoro HD, Wirya IGNW, Pudjiaji AH, Bisanto J, Zulkarnaen SZ, Uji diagnostik. Dalam: Sastroasmoro S, editor. Dasar-dasar Metodologi Penelitian Klinis. Edisi ke-3. Jakarta: Sagung Seto; 2010: 221-43
49. Altman DG, Bland JM. Diagnostic test 1: Sensitivity and specificity. 1994. BMJ 308 (6943): 1552

50. Akobeng AK. Understanding diagnostic test 2: likelihood ratios, pre-and post-test probabilities and their use in clinical practice. *Acta Paediatr.*2007; 96(4): 487-91
51. Fegan DF. Evaluation of Diagnostic Test: The Epidemiological Approach. Diunduh tanggal 9 April 2012 dari:
<http://www.fao.org/DOCREP/005/X4946E/x4946e0b.htm>
52. Mirsal Tawi. Wanita penjaja sek, pelanggannya dan HIV/AIDS.2009. Available at URL: <http://syehaceh.wordpress.com>. Diakses tanggal 3 November 2011
53. Paramanik JM, Potdar S, Kerkar S. Diagnosis of *Chlamydia trachomatis* infection. *J Clin Lab Anal* 2006; 20: 8-14
54. Krishnaputri S. Keandalan *Trichomonas rapid test* untuk Diagnosis Trikomoniasis Wanita Pekerja Seks Komersial Di Klinik Graha Sriwijaya Palembang. Tesis. Palembang: Bagian Ilmu Kesehatan Kulit dan Kelamin FKUNSRI, 2011
55. Philips AJ. Chlamydial infections. In: Nelson AL, Woodward JA, eds. A practical Guide for Primary Care Sexually Transmitted Diseases. 1st Totowa: Humana Press; 2007: 127-50
56. Sedyaningsih ER, Tanudyaya FK, Anartati AS, Purnamawati KA, Sutrisna A, Siswadi dkk. Laporan hasil penelitian prevalensi infeksi saluran reproduksi pada wanita penjaja seks di Palembang, Sumatera Selatan Indonesia, 2005. : 1-54
57. Debatista J, Martin D, Jamieson J, Crane K, Dalton I, Hall SR. et.al. Detection of *Chlamydia trachomatis* in an Australian high school student population. *Sex Transm infect* 2002; 78: 194-97