


SeroFIA™ Chlamydia

CE 0483



The SeroFIA™ Chlamydia product line consists of IgG, IgM and IgA specific immunofluorescence-based assays for the differential determination of Chlamydia: pneumoniae, trachomatis, and psittaci antibodies in human serum.

The user-friendly SeroFIA™ Chlamydia slides contain one row of seven antigen-coated wells for each of the three Chlamydia species.

- ▶ **Accurate**
Highly sensitive and specific
- ▶ **User freindly**
Simultaneously differentiates the three chlamydia species: *C. trachomatis*, *C. pneumoniae* and *C. psittaci*
Convenient design reduces error
Kit includes all reagents and controls
- ▶ **Reliable**
Excellent correlation with other commercial and in-house MIF assays

Ordering information

No. of tests/kit	Product name	Catalog number
105	SeroFIA™ IgG Chlamydia	511-01
105	SeroFIA™ IgM Chlamydia	512-01
105	SeroFIA™ IgA Chlamydia	513-01
105	SeroFIA™ <i>C. psittaci</i>	570-01
105	SeroFIA™ <i>C. trachomatis</i>	580-01
105	SeroFIA™ <i>C. pneumoniae</i>	590-01



SEROFIA™