

# CoproStrip™

for detection of *Giardia* or/and *Cryptosporidium*


**The newly designed CoproStrip™** Device tests are rapid chromatographic immunoassays for the qualitative detection of *Cryptosporidium* or/and *Giardia* antigens in human faecal specimens to aid in the diagnosis of *cryptosporidiosis* and *giardiasis*.

## Advantages:


- ♥ Newly designed Device- **Cassette format**
- ♥ Clear visual results: **Green/Red/Blue** colored test/control lines
- ♥ **Superior sensitivity and specificity**
- ♥ **Stool antigen extraction device** - Included
- ♥ **Results in just a few minutes**



## Test procedure:



## Interpretation of results:



## Test's performance:

Product name	Sensitivity	Specificity
CoproStrip™ Giardia	>99%	>99%
CoproStrip™ Cryptosporidium	>99%	>99%
CoproStrip™ Giardia/Cryptosporidium	98%	>99%

\* Reference: Microscopy

## Ordering information

Kit name	Tests/Kit	Catalog No.
CoproStrip™ Giardia	20	41217
CoproStrip™ Cryptosporidium	20	41218
CoproStrip™ Giardia/Cryptosporidium	20	41219